

## ACCESSIBILITY EQUIPMENT – updated July 2020



### Escape Lite 338-S attendant propelled wheelchair

A general use wheelchair which has a folding frame and back rest, which folds to a width of 28cm/11" with a standard 46cm/18" seat. There are attendant cable operated brakes and convenient swing away or detachable footplates. Maximum user weight 16 stone (101½ kg).




### Breezy self or attendant propelled general use wheelchair

This model has extra frame holes to allow for the height adjustable seat (maximum seat height is 51cm/20"). There are adjustable footplates and the seat size is 43x43cm (17X17"). Maximum user weight 19 ½ stone (125kg).



### Coopers self or attendant propelled shower wheelchair

This can be used as either a showering wheelchair or a commode chair (if used with a commode pan). It has a height adjustable frame and detachable or height adjustable footrests, also for extra comfort there is a padded seat which includes a removable section for easier personal cleaning. There are also padded armrests that can fold up out of the way. Seat height 54 - 63.5cm (21¼ - 25"). Seat dimensions 43x43cm (17X17").



### Wooden Transfer Board

Made from sturdy birch plywood this transfer board can be used to assist with transferring from beds, wheelchairs and cars without the need to be lifted. The board has tapered edges making it easier for the user to position themselves and the lacquered finish to the top gives a slippery surface to aid in sliding. Maximum user weight 30 stone (190kg).

### Suspended bath seat



This removable/repositionable seat provides a sitting platform within the bath to facilitate washing or can provide a staged decent into the bath. (Note: a user needs to be able to safely step into/out of the bath and position themselves. This therefore requires good upper body strength; Grab bars could be positioned by the bath to assist as required by individual guests).

### Cone furniture raisers



Raises a single bed\* or chair by 9cm/3 ½" and holds furniture legs and coasters in a firm position. This aids movement in getting in and out of furniture and reduces the strain on the body. Maximum user and load weight when 4 raisers are used is 35 stone (225kg). The raisers can also be used in conjunction with a powered hoist, which can be hired by prior arrangement.

\* Not suitable for a double bed configuration.



#### Extendable suction grab bars

The grab bars are easy to install or remove, just by flipping the tab. Two sizes are available 44cm (17") and 55cm (22") and both are extendable by up to 13cm (5"). They provide a safe means of getting in or out of the shower or bath and can withstand about double the pull pressure that a person could physically achieve.



#### Days wheeled adjustable height walking aid

Due to a lightweight aluminium frame and easy to manoeuvre wheels on the front, this walker can give confidence and stability. There are also slip resistant rubber tips on back and plastic handles for added safety.



#### Lightweight wheeled rollator with rest seat

This walker has 4 easily manoeuvrable wheels and a lightweight frame. It has an adjustable handle height for added comfort and there is a padded seat and backrest for taking breaks. Maximum user weight 18 ½ stone (117kg).



### Bed rails

These adjustable metal bed rails can be fitted to the side of the bed to prevent the risk of either falls or injuries, therefore giving the freedom to relax comfortably while in bed without the worry of slipping out. This can also aid in sitting up or changing of position without assistance. Height 46cm (18") Length 127cm (50").



### Bed Grab Rail

This sturdy bed grab rail ensures extra assistance when getting into or out of bed. With a height adjustable frame, the grab rail is suitable for any user. The two legs rest on the floor and the frame is firmly secured into place via three straps around the bed frame, ensuring minimal movement when in use.



### Vibrating alarm clock

This clock has several features to make getting up easier as it has not only a super-loud adjustable volume (up to 95db), but also an extra bright strobe alarm and a large, clear digital display. There are options to use as a shaker only, the shaker with sound, the sound and flash or the flash and shaker.



### Homecraft raised toilet seat with frame

The raised toilet seat is affixed to a sturdy tubular frame with moulded rubber slip resistant feet which provides extra security when using the toilet. Frame height is adjustable. Supports up to 25 stone (158kg).



#### Homecraft commode

Light and durable steel frame with adjustable height. The rear framework adds support for the back and arm rests provide extra comfort for the user. This offers a solution to those with limited mobility as it can be used anywhere. Maximum user weight 21 stone (136kg).



#### Raised toilet seat with frame

Height adjustable frame, which is designed to support the user in a high position for easy transfer on and off the seat. There are also rubber moulded feet to help prevent the frame from slipping. It has a wide 46cm (18") seat. Maximum user weight 30 stone (190 kg).



#### Prima Modular perching stool with arms and back

The tilt design seat is 5cm/2" lower at the front than the back, which reduces the effort required to sit/stand and offers an arm and back rest for added comfort. There is an adjustable seat height with a minimum height of 48cm/19" and a maximum height of 63cm/25" at the front. The user can half stand or perch, which in turn reduces the strain from standing up for long durations. The seat is 40cm/16" wide and the width between arms is 43cm/17". Maximum user weight 30 stone (190kg).



#### Coopers elbow crutch

Height adjustable elbow crutch with cuff around forearm and hand grip. The shape is designed to spread the pressure more evenly over the hand and therefore makes it more comfortable to use. Maximum user weight 28 stone (180kg).



#### Hand Bell

A convenient and loud means of attracting attention. It can be moved around from room to room and can therefore be placed either at the side of the bed, a chair in the sitting room or can be taken into the bathroom to provide extra confidence for personal independence.

**The following items are available to be hired from Solent Mobility  
at your request prior to arrival**

**for further details go to [Solent Mobility.co.uk](http://SolentMobility.co.uk)**



**Birdie hoist**

This electronically operated hoist is safe and easy to use, with its long reach and lots of leg space it enables the user to be rotated to 360 degrees in comfort and allows the user to be lifted from either a bed, a chair or the floor. Maximum user weight 23.5 stone (149kg).



**Profiling bed**

An electrically operated bed which offers personalised comfort and support to the user during sleep, activities and mealtimes. The optimised dimensions of the sections of the mattress platform are designed to reduce shear and friction, which ensures comfort whilst preventing migration down the bed. Maximum user weight 26 stone (165kg).



**Pearl mobility scooter\***

A lightweight and simple to control scooter that can be disassembled easily giving greater flexibility for exploring. A fully charged battery will allow 15 miles usage at 4mph. Padding throughout provides comfort whilst the swivel seat and flip-up armrests make accessing the scooter easier. Maximum user weight 21 ½stone (137kg).

\*Other styles of scooter may be available on asking. Scooter hire and shop mobility is also available at Petersfield, Winchester, Southampton, Portsmouth and in many other nearby towns.